Ready Go[™]

ReadyGo Diagnostics Ltd Newark Works, Bath Quays Bath, BA2 3DQ UK

info@readygotest.com www.readygotest.com

Snapshot

The complete solution – sample collection, preparation and molecular testing

- Small portable device, designed for consumer home-based use and professional, decentralised diagnostic testing
- Simple sample collection procedure The swab and collection buffers are integrated into one disposable for ease of use
- Collect identical amount of sample each time
- Lyses cells, viruses and bacteria to expose more starting material and increase test sensitivity
- Designed to collect a wide range of biological samples including saliva, buccal and nasal specimens

- Reduces sample debris and eliminates sample centrifugation steps
- No requirement to purify RNA and DNA for molecular tests
- Results within 20 minutes of sampling
- Capable of conducting singleplex and multiplex tests
- USB-C power connection
- Secure data transfer, facilitated by NFC
- Conveniently view results on the app
- Cloud based data storage
- No additional laboratory equipment required

Transforming home-based and decentralised diagnostic testing

